


## PIKA HEATEx 36

For the rapid DNA isolation from enriched samples

Cat. No. 2034-1-1

 **Warning!** Read the manual and the Safety Data Sheets before starting the analysis. Safety Data Sheets are available in the download area from [www.pika-weihenstephan.de](http://www.pika-weihenstephan.de). All handling steps should be performed under sterile conditions. Wear appropriate protective clothing and powder-free gloves. The use of filter tips is recommended. For testing of Food and Environmental samples only.

### Product description

The DNA Isolation Kit PIKA HeatEx 36 is intended for the extraction of DNA from bacteria and yeast cells in enriched samples and is optimized for Real time PCR. The microorganisms can either be isolated from solid surfaces i.e. colonies on membrane filters, or out of liquid samples. Prior isolation of single colonies is not necessary. Extracted DNA can directly be used for molecular analyses like PCR.

### Kit content

Materials supplied are sufficient for 36 isolations.

| Description                   | Amount      | Storage*                |
|-------------------------------|-------------|-------------------------|
| Washing buffer A (yellow cap) | 1 x 10.0 mL | 4°C or room temperature |
| Lysis buffer B (blue cap)     | 1 x 10.0 mL |                         |
| Dilution buffer A             | 1 x 2,0 mL  |                         |

\* Kit is shipped at ambient temperature

### Materials required but not supplied

| Material   |
|--|
| <b>Instruments and equipment</b>                   |
| Benchtop microcentrifuge for 1.5 mL reaction tubes |
| Reaction tube mixer (Vortexer)                     |
| Thermoincubator or water bath set to 80°C          |
| Pipettors  |
| <b>Consumables and reagents</b>                    |
| Powder-free gloves                                 |
| 1.5 mL reaction tubes, safe-lock, sterile          |
| Filter pipette tips                                |

## Procedural guidelines

1. Transfer the sample into a 1.5 mL reaction tube:
  - a) *Liquid samples:*
    - 50 µL of a turbid, bacterial sample (previously enriched sample or spoiled product)
    - 1.0 – 1.5 mL of a clear sample (even larger sample sizes can be used)
    - 50 – 200 µL of yeast slurries to reach a pellet size of app. 2 mm in diameter after centrifugation (see fig. 1)
  - b) *Colonies:* single colonies as well as different colonies can be processed together as one sample
    - Transfer 200 µL Washing buffer A and cell material into a 1.5 mL reaction tube. Skip steps 4. to 6.
2. Centrifuge for 3 min at 14,000 rpm (25,000 x g) or alternatively 10 min at 4,000 rpm (1,500 x g)
3. Control the pellet size. Pellet size should not exceed 2 mm in diameter (see fig. 1)  
If necessary, remove part of the pellet together with the liquid phase
4. Remove the liquid phase carefully and discard
5. Optional: Wash the pellet with 200 µL Washing buffer A, resuspend pellet and repeat steps 2 and 3
6. Add 200 µL of Lysis buffer B and resuspend by mixing briefly
7. Incubate sample at 80 °C ± 5 °C for 10 min in a thermoincubator or water bath
8. Centrifuge again as in step 2  
The pellet contains cell walls and other particles separated from the DNA
9. Transfer 100 µL of the liquid phase containing the DNA to a new 1.5 mL reaction tube and use the liquid phase for PCR.  
For long-term storage, freeze at -18 to -20 °C




Fig. 1: recommended pellet sizes:  
Left: max. pellet bacterial size  
Right: max. pellet size for samples containing yeast

Dilution buffer A may be used to dilute the DNA prior to further analysis.

PIKA Weihenstephan GmbH  
Raiffeisenstraße 31A  
85276 Pfaffenhofen  
GERMANY  
Phone +49(08441)879 48 30  
Fax +49(08441)879 48 31  
[www.pika-weihenstephan.de](http://www.pika-weihenstephan.de)  
[order@pika-weihenstephan.de](mailto:order@pika-weihenstephan.de)

**Notes: Use of product:** This product is to be used for research purposes only. **Property Rights:** The procedure used for the sample preparation as well as the kit are patent pending by PIKA Weihenstephan GmbH. For any commercial use of the kit or parts of this licensing from PIKA Weihenstephan GmbH is required. The use of our products may touch property rights of third parties. PIKA Weihenstephan GmbH does assume no responsibility for the lawful use of this kit; this responsibility lies expressly and solely at the user.